

Cardio-GeneSGKit®

Sistemas Genómicos has developed Cardio-GeneSGKit® series for the study of up to 235 genes associated with the diagnosis of several heterogeneous diseases, either syndromic or non-syndromic, whether autosomal dominant, autosomal recessive or X chromosome-linked. Thanks to the Cardio-GeneSGKit®, a broad range of pathologies are being studied such as familial cardiopathies (channelopathies, coronary heart diseases, cardiomyopathies), aortic pathologies and sudden death.

The Cardio-GeneSGKit® series boasts CE-IVD approval and contains reagents for processing twelve samples with the corresponding bio-computing analysis through the GeneSystems® platform. The workflow is optimised so as to carry out the libraries in nine hours and the technology used allows the combination of multiple samples and different GeneSGKits® in a single run on Illumina® sequencing platforms, ensuring an equitable balance of the data among the different samples.

HIGHLIGHTS

The most complete solution on the market

The Cardio-GeneSGKit® is cleared for CE-IVD and contains reagents to process twelve samples with the corresponding bioinformatics analysis and visualization by the GeneSystems® platform.

Fast workflow and low input requirements


The libraries are ready within nine hours of work starting with 50 ng of DNA.

Complete bioinformatic analysis and data visualization

The bioinformatics pipeline is executed in less than 48 hours, providing the user with a rich variant annotation including public, private and Sistemas Genómicos' reference databases. This comprehensive annotation reports SNPs and indels, which are visualized and easily filtered on GeneSystems®.

CNVs detection

Sensitive CNVs detection validated in prestigious public and private organisations.



**GeneSGKits® library
preparation**




**Sample sequencing on
Illumina®* platforms**



**Analysis and data interpretation on
GeneSystems®**

*HiSeq and MiSeq. Other sequencing platforms in validation

GeneSGKit® WORKFLOW



ABOUT US

Sistemas Genómicos is a pioneering company in the application of Next Generation Sequencing technology for genetic diagnosis. Accredited and certified by the major entities in quality assurance (ENAC, AENOR and CLIA), we offer solutions to accelerate and standardize clinical genetic testing, enabling clinicians to provide accurate diagnostics to their patients.

OUR MISSION

We strive to provide the best diagnostic solutions with the best technology and knowledge available to maximise patient benefit. Our integrated systems have helped laboratories for the last two decades to improve their diagnostics and research work, by developing a unique range of CE-IVD solutions.



DESIGN

The GeneSGKits® offer excellent performance from only 50 ng of DNA. The capture probes are based on ultra-long RNA baits (120 nucleotides) with high sensitivity and specificity for the detection of SNPs, indels and CNVs. Our probe design captures all the coding regions of each targeted gene, including splicing sequences, with a recommended average coverage of 200x.



PERFORMANCE

With GeneSGKits® clinicians and researchers can be confident in the accuracy and reproducibility of the sequencing data. Each sample is sequenced with a minimum coverage of 20x in over 98% of the regions of interest for relevant pathologies, guaranteed by strict validation protocols.



SUPPORT

For Sistemas Genómicos, customer is at the center of our work. This includes a personalised service to enable a seamless implementation of diagnostic next generation sequencing. We offer personalised training programs to ensure a quick acquisition of the skill set required to successfully and reliably harness the latest advances in sequencing technology to maximize patient benefit.

GeneSGKit®

GeneSystems® is a leading platform for the analysis of NGS data for diagnosis and genetic research which is registered as a health product with a CE Marking for diagnostic use in vitro for processing information from massive DNA sequencing to detect genetic variants and help to diagnose at the molecular level heterogeneous genetic diseases by means of directed resequencing. Developed in a private, secure and reliable Cloud Computing environment, it includes:



**Uploading and
downloading data**



**Complete access to all the data, from
raw data to annotated variants**



**Visualization and
prioritization of results**

It is developed by a multidisciplinary team of geneticists, bioinformatic and engineers, and validated by more than 200 companies, including hospitals, laboratories, research centers, universities, pharmaceutical companies. The platform has been audited independently and allows the detection of SNP and INDELS.

DATA ANALYSIS

GeneSystems® is designed to solve the different stages of the analytical data cycle of massive sequencing and provides solutions from the FASTQ files to the clinically-relevant data.

GeneSystems® captures annotations at gene, variant and transcript level and incorporates the most clinically relevant databases for the interpretation of variants, such as Ensemble, 1000Genomes, HGMD, dbSNP, RefSeq, Interpro, Prosite, Cosmic, HapMap, OMIM, ClinVar, ExAC, etc. Furthermore, GeneSystems® users have access to the proprietary variant database of Sistemas Genómicos. Within the platform there is also information on conservation scores, pathogenicity in silico predictors, flanking regions, among others. A one stop solution to find clinical meaning in your sequencing data.

GeneSystems® provides a highly sensitive CNVs detection validated in prestigious public and private organisations, which optimizes the MLPA experiments thus reducing costs and diagnostic turnaround times.

Hight data quality.

With the GeneSGKits® technology, clinicians and researchers will have the best tools for genetic diagnosis, thanks to the excellent quality of the sequencing data of the samples prepared with GeneSGKit®. Each sample will be sequenced with low variance in the coverage values with 98% of the regions sequenced at a minimum depth of 20x.

Technical data

GeneSGKit®	Genes	Size (MB)	% COVERAGE	
			12x	20x
Cardio-GeneSGKit®	235	1	99.4712	99.3306
Cardio-GeneSGKit® Arritmias	62	0.49	99.97	99.57
Cardio-GeneSGKit® MCP	90	0.75	100	99.91
Cardio-GeneSGKit® CVascular	54	0.35	99.78	99.41
Cardio-GeneSGKit® DMM	74	0.82	99.41	98.99

References

Cardio-GeneSGKit®: LV3891 (CE-IVD) / LV3892 (RUO)
 Cardio-GeneSGKit® Arritmias: LV3454 (CE-IVD) / LV2554 (RUO)
 Cardio-GeneSGKit® MCP: LV3460 (CE-IVD) / LV2874 (RUO)
 Cardio-GeneSGKit® CVascular: LV3455 (CE-IVD) / LV2558 (RUO)
 Cardio-GeneSGKit® DMM: LV3461 (CE-IVD) / LV2877 (RUO)

Genes description

Cardio-GeneSGKit®

A2ML1	CAV1	EFEMP2	GATA6	KCNE1	MIB1	OBSCN
ABCC6	CAV3	EHMT1	GATAD1	KCNE2	MMP21	PDLIM3
ABCC9	CBL	EIF2AK4	GBE1	KCNE3	MYBPC3	PKP2
ACTA2	CBS	ELN	GDF1	KCNE4	MYH11	PKP4
ACTC1	CFAP53	EMD	GDF2	KCNE5	MYH6	PLEC
ACTN2	CFC1	ENG	GJA1	KCNH2	MYH7	PLN
ACVR2B	CHD7	ENPP1	GJA5	KCNJ2	MYL2	PLOD1
ACVRL1	CHST14	EVC	GLA	KCNJ5	MYL3	PRDM16
AGL	CITED2	EYA4	GLMN	KCNJ8	MYLK	PRKAG2
AKAP9	COL1A1	FBLN5	GNAI2	KCNK17	MYLK2	PRKG1
ALG10	COL1A2	FBN1	GPD1L	KCNK3	MYOT	PROC
ANK2	COL3A1	FBN2	GUCY1A3	KCNQ1	MYOZ2	PROS1
ANKRD1	COL4A1	FHL1	GYG1	KRAS	MYPN	PSEN1
BAG3	COL5A1	FHL2	GYS1	KRIT1	NEBL	PSEN2
BMPR2	COL5A2	FHOD3	HAND2	LAMA4	NEXN	PTPN11
BRAF	CRELD1	FKBP14	HCN4	LAMP2	NF1	RAF1
CACNA1B	CRYAB	FKRP	HFE	LDB3	NKX2-3	RANGRF
CACNA1C	CSRP3	FKTN	HRAS	LDLR	NKX2-5	RASA1
CACNA1D	CTNNA3	FLNA	HTRA1	LMNA	NKX2-6	RASA2
CACNA2D1	DES	FLNC	ILK	LOX	NODAL	RBM10
CACNB2	DMD	FOXC1	JAG1	LRP6	NOTCH1	RBM20
CALM1	DPP6	FOXF1	JPH2	LZTR1	NOTCH2	RIT1
CALM2	DSC2	FXN	JUP	MAP2K1	NOTCH3	RNF213
CALM3	DSG2	G6PC3	KCNA5	MAP2K2	NPPA	RPSA
CALR3	DSP	GAA	KCND2	MEF2A	NRAS	RRAS
CASQ2	DTNA	GATA4	KCND3	MFAP5	NUP155	RYR1

RYR2	SHOC2	SNTA1	TBX1	TGFB3	TNNI3K	VCL
SCN1B	SLC25A4	SOS1	TBX20	TGFBR1	TNNT2	XK
SCN2B	SLC2A10	SOS2	TBX4	TGFBR2	TNXB	ZFPM2
SCN3B	SLMAP	SPRED1	TBX5	TLL1	TPM1	ZIC3
SCN4B	SMAD3	SYNE1	TCAP	TMEM43	TRDN	ZNF469
SCN5A	SMAD4	SYNE2	TEK	TMPO	TRPM4	
SDHA	SMAD6	TAB2	TFAP2B	TNNC1	TTN	
SGCD	SMAD9	TAZ	TGFB2	TNNI3	TTR	

Cardio-GeneSGKit® Arritmias is designed to diagnose the following pathologies: Brugada syndrome, long QT syndrome, short QT syndrome, atrial fibrillation, catecholaminergic polymorphic ventricular tachycardia and arrhythmogenic right ventricular dysplasia.

Genes description

Cardio-GeneSGKit® Arritmias

ABCC9	CASQ2	GJA1	KCNE1L	KCNJ5	PNN	SCN5A
AKAP9	CAV3	GJA5	KCNE2	KCNJ8	RANGRF	SCNN1B
ANK2	CRYAB	GPD1L	KCNE3	KCNQ1	RPSA	SCNN1G
BAG3	DES	HCN1	KCNE4	KCNQ2	RYR2	SNTA1
CACNA1B	DPP6	HCN4	KCNH2	LDB3	SCN10A	TGFB3
CACNA1C	DSC2	JUP	KCNJ11	MYOT	SCN1B	TMEM43
CACNA1D	DSG2	KCNA5	KCNJ12	NPPA	SCN2B	TRPM4
CACNA2D1	DSP	KCND3	KCNJ2	PKP2	SCN3B	TTN
CACNB2	FLNC	KCNE1	KCNJ3	PKP4	SCN4B	

Cardio-GeneSGKit® MCP is designed to diagnose the following pathologies: hypertrophic, dilated and restrictive cardiomyopathies, ventricular non-compaction, cardiomyopathies secondary to muscular dystrophies, arrhythmogenic ventricular dysplasia, RASopathy, Noonan syndrome, Legius syndrome and cardiofacial syndrome.

Genes description

Cardio-GeneSGKit® MCP

ABCC9	CBL	FHL2	LDB3	NEBL	RAF1	TAZ
ACTC1	CRYAB	FKTN	LMNA	NEXN	RBM20	TCAP
ACTN2	CSRP3	FLNC	MAP2K1	NRAS	RPSA	TGFB3
ADRB1	CTF1	FXN	MAP2K2	PDLIM3	RYR2	TMEM43
ADRB2	DES	GAA	MYBPC3	PKP2	SCN5A	TMPO
ADRB3	DMD	GLA	MYH6	PKP4	SDHA	TNNC1
AGL	DSC2	HRAS	MYH7	PLEC	SGCD	TNNI3
ANK2	DSG2	ILK	MYL2	PLN	SHOC2	TNNT2
ANKRD1	DSP	JPH2	MYL3	PNN	SLC25A4	TPM1
BAG3	DTNA	JUP	MYLK2	PRKAG2	SOS1	TTN
BRAF	EMD	KRAS	MYOT	PSEN1	SPRED1	TTR
CALR3	EYA4	LAMA4	MYOZ2	PSEN2	SYNE1	VCL
CAV3	FHL1	LAMP2	MYPN	PTPN11	SYNE2	

Cardio-GeneSGKit® CVascular is designed to diagnose the following pathologies: aortic aneurysm, hypercholesterolaemia, monogenic coronary diseases and congenital cardiopathies.

Genes description

Cardio-GeneSGKit® CVascular

ACTC1	APOE	ELN	GDF1	MEF2A	PCSK9	TGFB1
ABCA1	ANGPTL3	ENPP1	GJA1	MYH11	PLD3	TGFB2
ABCC6	BMPR2	FBLN5	GLMN	MYH6	RASA1	TGFB3
ABCG5	CFC1	FBN1	JAG1	MYLK	RBM10	TLL1
ABCG8	COL3A1	FBN2	LDLR	NKX2-3	SLC2A10	ZNF469
ACTA2	COL4A1	G6PC3	LDLRAP1	NKX2-5	SMAD3	ZFPM2
APOA1	CYP7A1	GATA4	LPL	NKX2-6	TBX1	
APOB	EFEMP2	GATA6	LRP6	NOTCH1	TBX20	

Cardio-GeneSGKit® DMM is designed to diagnose different types of muscular dystrophies and myopathies that follow autosomal dominant and X-linked hereditary patterns.

Genes description

Cardio-GeneSGKit® DMM

ACTA1	CLCN1	DOK7	ISCU	MYOT	RYR1	TPM3
AGRN	CNTN1	DYSF	ITGA7	NEB	SEPN1	TRIM32
ANO5	COL6A1	EMD	KBTBD13	OPA1	SGCA	TRPV4
ATP2A1	COL6A2	FHL1	LAMA2	PABPN1	SGCB	TTN
ATP7A	COL6A3	FKRP	LARGE	PLEC	SGCG	UBA1
BAG3	COL9A3	FKTN	LDB3	PLEKHG5	SYNE1	VAPB
BIN1	CRYAB	FLNC	LMNA	POMGNT1	SYNE2	VCP
CAPN3	DAG1	GNE	MATR3	POMT1	TAZ	VRK1
CAV3	DES	GSN	MTM1	POMT2	TCAP	
CFL2	DMD	HSPG2	MYH2	PSMB8	TNNT1	
CHKB	DNM2	IGHMBP2	MYH7	PTRF	TPM2	

