

LMA-GeneSGKit[®]

Sistemas Genómicos has developed LMA-GeneSGKit[®] for the study of 62 genes involve in myeloid malignances. LMA-GeneSGKit[®] allows the diagnosis and monitoring of different pathologies such as Acute myeloid leukemia (AML), Myelodysplastic syndrome (MDS), and Chronic Myeloproliferative disorders.

The LMA-GeneSGKits[®] are cleared for CE-IVD and contain reagents for processing twelve samples with the corresponding bioinformatic analysis on the GeneSystems[®] platform. Developed in collaboration with experts myeloid malignances, the LMA-GeneSGKit[®] is designed for the study of 62 genes associated with Acute myeloid leukemia (AML), Myelodysplastic syndrome (MDS), and Chronic Myeloproliferative disorders. The workflow is optimised for library preparation in nine hours and the technology used allows the combination of multiple samples prepared with different GeneSGKits[®] in a single run on Illumina[®] sequencing platforms

HIGHLIGHTS

The most complete solution on the market

The LMA-GeneSGKit[®] is cleared for CE-IVD and contains reagents to process twelve samples with the corresponding bioinformatics analysis and visualization by the GeneSystems[®] platform.

Fast workflow and low input requirements


The libraries are ready within nine hours of work starting with 50 ng of DNA.

Complete bioinformatic analysis and data visualization

The bioinformatics pipeline is executed in less than 48 hours, providing the user with a rich variant annotation including public, private and Sistemas Genómicos' reference databases. This comprehensive annotation reports SNPs and indels, which are visualized and easily filtered on GeneSystems[®].

High-accuracy and high-sensitivity analysis

Limit of detection down to 10% mutant allele frequency with 100x minimum coverage of each region.



GeneSGKits[®] library preparation




Sample sequencing on Illumina[®]* platforms



Analysis and data interpretation on GeneSystems[®]

*HiSeq and MiSeq. Other sequencing platforms in validation

GeneSGKit® WORKFLOW



ABOUT US

Sistemas Genómicos is a pioneering company in the application of Next Generation Sequencing technology for genetic diagnosis. Accredited and certified by the major entities in quality assurance (ENAC, AENOR and CLIA), we offer solutions to accelerate and standardize clinical genetic testing, enabling clinicians to provide accurate diagnostics to their patients.

OUR MISSION

We strive to provide the best diagnostic solutions with the best technology and knowledge available to maximise patient benefit. Our integrated systems have helped laboratories for the last two decades to improve their diagnostics and research work, by developing a unique range of CE-IVD solutions.



DESIGN

The GeneSGKits® offer excellent performance from only 50 ng of DNA. The capture probes are based on ultra-long RNA baits (120 nucleotides) with high sensitivity and specificity for the detection of SNPs and indels. Our probe design captures all the coding regions of 42 genes including splicing sequences plus exonic hotspots of an additional 20 genes.



PERFORMANCE

With GeneSGKits® clinicians and researchers can be confident in the accuracy and reproducibility of the sequencing data. LMA-GeneSGKit® provides a 100x minimum coverage of each region at 300x in over 93% of the region of interest guaranteed by strict validation protocols.



SUPPORT

For Sistemas Genómicos, customer is at the center of our work. This includes a personalised service to enable a seamless implementation of diagnostic next generation sequencing. We offer personalised training programs to ensure a quick acquisition of the skill set required to successfully and reliably harness the latest advances in sequencing technology to maximize patient benefit.

 GeneSGKit®



GeneSystems®

SISTEMAS GENÓMICOS

GeneSystems® is a leading platform for the analysis of NGS data for diagnosis and genetic research which is registered as a health product with a CE Marking for diagnostic use in vitro for processing information from massive DNA sequencing to detect genetic variants and help to diagnose at the molecular level heterogeneous genetic diseases by means of directed resequencing. Developed in a private, secure and reliable Cloud Computing environment, it includes:



**Uploading and
downloading data**



**Complete access to all the data, from
raw data to annotated variants**



**Visualization and
prioritization of results**

It is developed by a multidisciplinary team of geneticists, bioinformatic and engineers, and validated by more than 200 companies, including hospitals, laboratories, research centers, universities, pharmaceutical companies. The platform has been audited independently and allows the detection of SNP and INDELS.

DATA ANALYSIS

GeneSystems® is designed to solve the different stages of the analytical data cycle of massive sequencing and provides solutions from the FASTQ files to the clinically-relevant data.

GeneSystems® captures annotations at gene, variant and transcript level and incorporates the most clinically relevant databases for the interpretation of variants, such as Ensemble, 1000Genomes, HGMD, dbSNP, RefSeq, Interpro, Prosite, Cosmic, HapMap, OMIM, ClinVar, ExAC, etc. Furthermore, GeneSystems® users have access to the proprietary variant database of Sistemas Genómicos. Within the platform there is also information on conservation scores, pathogenicity in silico predictors, flanking regions, among others. A one stop solution to find clinical meaning in your sequencing data.

GeneSystems® provides a highly sensitive CNVs detection validated in prestigious public and private organisations, which optimizes the MLPA experiments thus reducing costs and diagnostic turnaround times.

Hight data quality.

With the GeneSGKits® technology, clinicians and researchers will have the best tools for genetic diagnosis, thanks to the excellent quality of the sequencing data of the samples prepared with GeneSGKit®. Each sample will be sequenced with low variance in the coverage values with 98% of the regions sequenced at a minimum depth of 20x.

Technical data

			% COVERAGE
GeneSGKit [®]	Genes	Size (MB)	100x
LMA-GeneSGKit [®]	62	0,49	93,40*

*: % Coverage at 300x mean coverage

References

LMA-GeneSGKit[®] Ref.: LV3570 (Kit CE-IVD) / LV3569 (Kit RUO)

Genes description

LMA-GeneSGKit[®]

AKAP13	CREBBP	GNB1	MAZ	POU4F1	SHOC2	WT1 *
ASXL1 *	CSF3R*	HTT	MLL *	PPM1D	SMC3	ZRSR2
BCOR	CUX1	IDH1 *	MPL *	PTPN11 *	SOS1	
BCORL1	DNAH9	IDH2 *	MYD88	RAF1	SRSF2	
BRAF	DNMT3A	JAK2 *	NF1	RIT1	STAG1	
CA9	ELF4	KIT*	NOTCH1	RUNX1 *	STAG2	
CALR*	EZH2	KRAS*	NPM1 *	SETBP1	TET2	
CBL*	FLT3*	LTA4H	NRAS*	SETD2	TLR9	
CEBPA	FMN2	MAP1B	PHF6	SETDB1	TP53	
CNOT3	GNAS	MAP2K1	POU2F2	SF3B1 *	U2AF1 *	

*: Analysis of exonic hotspots