



LisosomalPeroxisomal-GeneSGKit®

Sistemas Genómicos has developed LisosomalPeroxisomal-GeneSGKit® with **145 genes whose mutations are associated with several lysosomal and peroxisomal disorders**. LisosomalPeroxisomal-GeneSGKit® allows the diagnosis of different pathologies such as lysosomal disorders (sphingolipidosis, mucopolysaccharidosis, mucopolipidosis and glycoproteinosis, lysosomal síndromes, lysosomal disorders with skeletal dysplasias, lysosomal disorders with metabolopathies, phospholipid disorders) and peroxisomal disorders (peroxisomal biogenesis disorder, defects of α and β -oxidation of fatty acids, other peroxisomal defects, Rhizomelic chondrodysplasia punctata, congenital disorders of glycation types I and II).

The workflow is optimised to prepare the libraries in nine hours and the technology used allows the combination of multiple samples and different GeneSGKits® in a single run on Illumina® sequencing platforms, ensuring an equitable balance of data from the different samples.

HIGHLIGHTS

The most complete solution on the market

LisosomalPeroxisomal-GeneSGKit® is cleared for CE-IVD and contains reagents to process twelve samples with the corresponding bioinformatics analysis and visualization by the GeneSystems® platform.

Fast workflow and low input requirements

The libraries are ready within nine hours of work starting with 50 ng of DNA.

Complete bioinformatic analysis and data visualization

The bioinformatics pipeline is executed in less than 48 hours, providing the user with a rich variant annotation including public, private and Sistemas Genómicos' reference databases. This comprehensive annotation reports SNPs and indels, which are visualized and easily filtered on GeneSystems®.

CNVs detection

Sensitive CNVs detection validated in prestigious public and private organisations.

1



GeneSGKits® library preparation

2



Sample sequencing on Illumina®* platforms


3



Analysis and data interpretation on GeneSystems®

*HiSeq and MiSeq. Other sequencing platforms in validation

GeneSGKit® WORKFLOW



ABOUT US

Sistemas Genómicos is a pioneering company in the application of Next Generation Sequencing technology for genetic diagnosis. Accredited and certified by the major entities in quality assurance (ENAC, AENOR and CLIA), we offer solutions to accelerate and standardize clinical genetic testing, enabling clinicians to provide accurate diagnostics to their patients.

OUR MISSION

We strive to provide the best diagnostic solutions with the best technology and knowledge available to maximise patient benefit. Our integrated systems have helped laboratories for the last two decades to improve their diagnostics and research work, by developing a unique range of CE-IVD solutions.



DESIGN

The GeneSGKits® offer excellent performance from only 50 ng of DNA. The capture probes are based on ultra-long RNA baits (120 nucleotides) with high sensitivity and specificity for the detection of SNPs, indels and CNVs. Our probe design captures all the coding regions of each targeted gene, including splicing sequences, with a recommended average coverage of 200x.



PERFORMANCE

With GeneSGKits® clinicians and researchers can be confident in the accuracy and reproducibility of the sequencing data. Each sample is sequenced with a minimum coverage of 20x in over 98% of the regions of interest for relevant pathologies, guaranteed by strict validation protocols.



SUPPORT

For Sistemas Genómicos, customer is at the center of our work. This includes a personalised service to enable a seamless implementation of diagnostic next generation sequencing. We offer personalised training programs to ensure a quick acquisition of the skill set required to successfully and reliably harness the latest advances in sequencing technology to maximize patient benefit.

 GeneSGKit®



GeneSystems®

SISTEMAS GENÓMICOS

GeneSystems® is a leading platform for the analysis of NGS data for diagnosis and genetic research which is registered as a health product with a CE Marking for diagnostic use in vitro for processing information from massive DNA sequencing to detect genetic variants and help to diagnose at the molecular level heterogeneous genetic diseases by means of directed resequencing. Developed in a private, secure and reliable Cloud Computing environment, it includes:



**Uploading and
downloading data**



**Complete access to all the data, from
raw data to annotated variants**



**Visualization and
prioritization of results**

It is developed by a multidisciplinary team of geneticists, bioinformatic and engineers, and validated by more than 200 companies, including hospitals, laboratories, research centers, universities, pharmaceutical companies. The platform has been audited independently and allows the detection of SNP and INDELS.

DATA ANALYSIS

GeneSystems® is designed to solve the different stages of the analytical data cycle of massive sequencing and provides solutions from the FASTQ files to the clinically-relevant data.

GeneSystems® captures annotations at gene, variant and transcript level and incorporates the most clinically relevant databases for the interpretation of variants, such as Ensemble, 1000Genomes, HGMD, dbSNP, RefSeq, Interpro, Prosite, Cosmic, HapMap, OMIM, ClinVar, ExAC, etc. Furthermore, GeneSystems® users have access to the proprietary variant database of Sistemas Genómicos. Within the platform there is also information on conservation scores, pathogenicity in silico predictors, flanking regions, among others. A one stop solution to find clinical meaning in your sequencing data.

GeneSystems® provides a highly sensitive CNVs detection validated in prestigious public and private organisations, which optimizes the MLPA experiments thus reducing costs and diagnostic turnaround times.

Hight data quality.

With the GeneSGKits® technology, clinicians and researchers will have the best tools for genetic diagnosis, thanks to the excellent quality of the sequencing data of the samples prepared with GeneSGKit®. Each sample will be sequenced with low variance in the coverage values with 98% of the regions sequenced at a minimum depth of 20x.

Lisosomal Peroxisomal-GeneSGKit® CATALOG

Technical data

GeneSGKit®	Genes	Size (MB)	% COVERAGE	
			12x	20x
Lisosomal Peroxisomal-GeneSGKit®	145	0.49	99.95	99.90

References

Lisosomal Peroxisomal-GeneSGKit® Ref.: LV3376 (Kit CE-IVD) / LV3377 (Kit RUO)

Genes description

Lisosomal Peroxisomal-GeneSGKit®

ABCD1	ASAH1	DNAJC5	GNPTG	MLYCD	PEX3	SPTLC1
ABHD12	ATP13A2	DNM1L	GNS	MOGS	PEX5	SPTLC2
ABHD5	ATP6V0A2	DPAGT1	GPC3	MPDU1	PEX6	SRD5A3
ACOX1	B4GALT1	DPM1	GUSB	MPI	PEX7	ST3GAL5
ACSL4	CAT	DPM3	HEXA	MVK	PHYH	SUMF1
ADAMTSL2	CHKB	DYM	HEXB	NAGA	PLA2G6	TAZ
AGA	CLN3	ELOVL4	HGSNAT	NAGLU	PMM2	TCF4
AGK	CLN5	FA2H	HSD17B4	NEU1	PNPLA6	TPP1
AGPS	CLN6	FUCA1	HYAL1	NPC1	PPT1	TRIM37
AGXT	CLN8	GAA	IDS	OFD1	PSAP	TUSC3
ALDH3A2	COG1	GALC	IDUA	PANK2	RAB33B	XDH
ALG1	COG7	GALNS	KDM6A	PEX1	RAI1	LAMP1
ALG12	COG8	GAMT	KMT2D	PEX10	RFT1	GBA2
ALG2	COL11A2	GATM	LAMP2	PEX11B	SCP2	GBA3
ALG3	COL2A1	GBA	LIPA	PEX12	SGSH	SCARB2
ALG6	CTNS	GLA	LPIN1	PEX13	SLC17A5	LYST
ALG8	CTSA	GLB1	MAN2B1	PEX14	SLC35A1	MYO5A
ALG9	CTSC	GM2A	MANBA	PEX16	SLC35C1	RAB27A
AMACR	CTSD	GNE	MCOLN1	PEX19	SLC6A8	MLPH
ARSA	CTSK	GNPAT	MFSD8	PEX2	SMPD1	
ARSB	DHCR7	GNPTAB	MGAT2	PEX26	SOD1	

